


Relax: ...ergonomic and functional

The stool matches with ComfortLine, UniversallLine and BasicLine chairs. The minimum height is 50 cm, maximum height is 70 cm. The seat can be tilted to the front. Five wheels with 5 cm diameter ensure high stability. Upholstery material and colours of applications are identical to those of the Bionic chairs.



Bi-colour design allows a colour combination fitting to the chair.

Standard colour:


Chrom grey 

Special colours for
bi-colour application:

Feuer 

Atoll 

Lima 

Amazonas 

Technical Data

Dimensions

Maximum height	70 cm
Minimum height	50 cm
Weight	7,5 kg
Wheels	Ø 5 cm



Multifunctional stool



Bionic Medizintechnik GmbH ... a JMS-corporation · Max-Planck-Straße 21
D-61381 Friedrichsdorf/Germany · Phone +49 (0)6172-7576-0
Fax +49 (0)6172-7576-10 · info@bionic-jms.com · www.bionic-jms.com

Right for technical modifications reserved
B-MOBRLEZ2





Relax

Multifunctional stool



The stool can be moved up and down and the seat can be tilted to the front.



Stool in maximum height



Stool in minimum height